

# HER-SPEC<sup>®</sup>

Single-Use **Click-Free** Vaginal Speculum with Integrated LED Light Source



Introducing **HER-SPEC**, a single-use, **CLICK-FREE** vaginal speculum with integrated LED light source. Developed by obp using decades of user feedback and experience, **HER-SPEC** improves the patient experience by eliminating the clicking sound and discomfort associated with the ratchet mechanism in disposable vaginal specula.

Ideal for clinics, emergency rooms, urgent care centers, and physician offices; **HER-SPEC** is designed for use on women undergoing a procedure requiring vaginal access and exposure.

Give your patients the best experience.

**\*HER-SPEC will be available for purchase in May of 2019**

## TO ORDER

EMAIL: [SALES@OBPMED.COM](mailto:SALES@OBPMED.COM)

CALL: 1-978-291-6853

WEB: [OBPMEDICAL.COM](http://OBPMEDICAL.COM)

- Item # C020500 (Small)
- Item # C020510 (Medium)
- Item # C020520 (Large)
- Item # C020600 (Small - Sterile)
- Item # C020610 (Medium - Sterile)
- Item # C020620 (Large - Sterile)

## ILLUMINATE THE EVERYDAY

**obp** is the leading global developer of single-use, self-contained, illuminating medical devices. We integrate light with device — making them one — creating effective and simplified single-use solutions. Our innovative products are used in more than 10,000 physician offices, surgery centers and hospitals throughout the U.S., as well as healthcare facilities worldwide. **At obp, we illuminate the everyday.**

# CLICK-FREE FOR HER COMFORT & EXPERIENCE



FOR MORE INFORMATION VISIT

[WWW.OBPMEDICAL.COM/CLICKFREE](http://WWW.OBPMEDICAL.COM/CLICKFREE)



360 Merrimack Street Lawrence, MA 01843  
978.291.6853 | [www.obpmedical.com](http://www.obpmedical.com)

REV 040219

## PRODUCT FEATURES

- Disposable, single-use device reduces risk of cross-contamination
- Each HER-SPEC comes ready for use with an integrated ultra-bright cordless LED light source
- Single-use, battery-operated LED light source can produce peak light for 30+ minutes
- Click-free engagement improves patient experience by eliminating the clicking sound and discomfort associated with common disposable vaginal specula ratchet mechanism
- Smooth locking feature allows clinicians to set blades to desired open position by pushing center tab to automatically lock blades in place
- Easily release blades by pushing tabs located to the right or left of center tab
- Ergonomically designed handle with diamond grip for added comfort and maneuverability during exams
- 105-degree angled handle creates extra space between clinician's hand and patients' body and/or exam table
- Battery removal feature allows users the option to easily eject batteries without contact after use
- No assembly required and no additional parts needed
- Also available, **ER-SPEC** single-use vaginal speculum with integrated LED light source and **OFFICESPEC** single-use side-opening vaginal speculum with integrated LED light source

Contract availability coming soon