

ONETRAC[®]

Single-Use Cordless Retractor with Integrated LED Light Source



obp[™]
simple, bright ideas.

Visualization, retraction and hemostasis during a surgical procedure is critical in reducing complications and infections. Avoid the costly and complicated reprocessing requirements of fragile fiber-optic cables. With the award winning **ONETRAC**, cordless access and illumination of the entire pocket is easily and safely achieved.

ONETRAC is an immediate and cost-effective solution for surgical procedures requiring access, illumination and exposure of a surgical tissue pocket. Illumination in the depths of a surgical tissue pocket or cavity can be met within seconds – no assembly, additional parts/components, or fiber-optic cables required.

TO ORDER

EMAIL: SALES@OBPMED.COM

CALL: 1-978-291-6853

WEB: OBPMEDICAL.COM

- Item # C100190 (175mm x 30mm)
- Item # C100110 (135mm x 30mm)
- Item # C100165 (135mm x 20mm)
- Item # C100180 (105mm x 60mm)
- Item # C100140 (90mm x 22mm)
- Item # C100170 (55mm x 8mm)
- Item # C100150 (40mm x 20mm)
- Sold in boxes of 5 EA/BX
- Each unit comes packed sterile

ILLUMINATE THE EVERYDAY

obp is the leading global developer of single-use, self-contained, illuminating medical devices. We integrate light with device — making them one — creating effective and simplified single-use solutions. Our innovative products are used in more than 13,000 physician offices, surgery centers and hospitals throughout the U.S., as well as healthcare facilities worldwide. **At obp, we illuminate the everyday.**

CUT THE CORD

New 175mm x 30mm

The first multi-LED ONETRAC



New 55mm x 8mm Aufricht

The first single-use Aufricht with integrated LED light source



PRODUCT FEATURES

- Disposable, single-use device reduces risk of cross-contamination
- On/off switch to preserve and maximize duration of peak illumination strength
- Integrated cordless LED light source produces no damaging heat and will provide peak power for 100 + minutes
- Integrated smoke evacuation channel
- Revolutionary reinforced nylon material provides exceptional strength and rigidity
- Angled tip with blunt distal end to maximize tissue adherence
- Lightweight for added comfort and ease of use

TIRED OF EXPENSIVE LIGHTED RETRACTOR SYSTEMS?

Today's lighted retractor systems cost thousands of dollars with hundreds more in ongoing disposable and/or reprocessing costs.

Each **ONETRAC** comes ready for use with everything needed to successfully retract soft tissue and illuminate a surgical pocket or cavity all at a fraction of the cost of other lighted retractor systems. No assembly, additional parts/components, or fiber-optic cables are required.

GPO CONTRACTS

VIZIENT CONTRACT MS3250

PREMIER CONTRACT PP-OR-1449

HEALTHTRUST CONTRACT 39624

INTALERE CONTRACT VH11348

ROI CONTRACT RMS02360

FIRSTCHOICE CONTRACT FC2370



360 Merrimack Street Lawrence, MA 01843
978.291.6853 | www.obpmedical.com

REV 010820